

Protocol: Tn5-on-beads DNA tagmentation

Version: v1.3, June 2020

By: Marek Kucka, Chan Lab, Friedrich Miescher Laboratory, Max Planck Campus Tübingen

Contact: mkucka@tuebingen.mpg.de; frank.chan@tue.mpg.de; <http://fml.tuebingen.mpg.de/chan-group.html>

Coupling of assembled Tn5 with magnetic beads.

Beads: 5ml of NEB, S1421, Hydrophilic Streptavidin Magnetic Beads

- Move 5 ml (or less/more) of NEB-beads from the bottle into a 50 ml falcon tube (shake the beads well before to re-suspend). Use a 15 ml falcon if using 2ml or less of beads, then put on magnetic stand and remove the bead-storage buffer.
- Add 20 ml of streptavidin binding buffer, mix to re-suspend the beads, put on magnet and remove the buffer.
- Add 30 ml of streptavidin binding buffer and the calculated amount of assembled Tn5.

Depending on the strength of the current batch of Tn5, e.g. I used 0.45 μ l Tn5 per 5 μ l NEB-beads => 450 μ l Tn5 in total for 5ml of NEB-beads in 30 ml of streptavidin binding buffer).

- Close the tube and immediately mix by inverting and let rotate the tube at 10 rpm at RT for 30 minutes.
- Put on magnet for 2-4 min and decant/remove the buffer.
- Add 30 ml of Dialysis buffer and rotate at 10 rpm at RT for 5 min.
- Put on magnet and remove buffer.
- Gently re-suspend the Tn5-on-beads in a desired amount of Dialysis buffer. I usually store the Tn5-on-beads in **4-times the original volume of NEB-bead** => 20 ml of dialysis buffer for every 5 ml of beads.
 - ➔ 10 μ l of Tn5-on-beads in dialysis buffer equals to 2.5 ul of original NEB-beads
 - ➔ 20 μ l equals to 5 ul of original NEB-beads.

Notes: For an easier resuspension of the beads, store the beads in a small bottle (e.g. glass 50 or 100 ml bottle), not a falcon tube. Store at 4C.

5 ml of original NEB beads => enough for 1000 samples when using 20ul of Tn5-on-beads per sample, or 2000 samples (if using 10 ul Tn5-on-beads per sample).

5ul of Tn5-On-beads cuts (tagments) maximum of approximately 8 ng of your input DNA (the rest will be discarded during washes).

Biotinylated-oligos used to assemble Tn5:

Biot-25xT-Tn5ME-A /5BiosG/TTCGTCGGCAGCGTCAGATGTGTATAAGAGACAG '3
Tn5ME-rev /5Phos/CTGTCTCTTATACACATCT '3

Biot-24xT-Tn5ME-B /5BiosG/TTGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG '3
Tn5ME-rev /5Phos/CTGTCTCTTATACACATCT '3

DNA Tagmentation with Tn5-on-beads

Transfer 20 µl or 10 µl of the Stock Tn5-on-beads per sample (20 µl = 5 µl of NEB beads; 10 µl equals to 2.5 µl of original NEB beads) into a 1.5 ml Eppendorf tube on magnet

- for a 96 samples experiment, you need 1 ml or 2 ml of Tn5-on-beads, then exchange the storage buffer for **1 ml of WASH buffer**. Keep on ice.

Note: Pre-warm the plate shaker and PCR-thermocycler to 55 °C.

Pipette 20 µl of tagmentation buffer per sample in a 96-well plate, then add DNA (1-4µl) and only then add 10 ul of Tn5-on-beads in Wash buffer:

1st: pipette 20 µl of tagmentation buffer in each well, per sample: **14 µl H2O**
 6 µl 5xTAPS-DMF-MgCl2 buffer

2nd: pipette 1-4 µl of gDNA to each well (but at least 5 ng per sample)

3rd: pipette 10 µl of Tn5-on-beads (in WASH buffer)

 Total = 30-34 µl

Seal and mix the plate in a plate shaker (>1500rpm), 30 sec, then incubate in PCR-thermocycler for 10 minutes at 55°C, then hold at 10°C.

Tip: put the plate on magnetic stand for 2 seconds and then straight into the plate-shaker if the beads have settled already and are difficult to re-suspend just by shaking in plate-shaker.

Note: shake the plate with at least 1500 rpm for 30 second (or by pipetting, or inverting the strip tubes if fewer samples) until the beads are evenly re-suspended.

OPTIONAL: mix in a plate shaker again after 5 minutes incubation to re-suspend the beads, and then continue the incubation at 55°C.

Stripping of the Tn5 from the DNA

Add 20 µl of "0.6% SDS- Wash" buffer to each sample.

Seal the plate and mix in the plate shaker or by inverting the tubes to re-suspend the beads.

Incubate at **55°C for 10 minutes**.

Put plate on magnet, remove supernatant.

Add 150 µl WASH buffer to the beads, incubate at RT for 1 minute, and then remove the wash buffer

Add 100 µl WASH buffer.

Keep in wash buffer until ready to pipet Q5 PCR-mix.

PCR amplification and adapter extension.

- Note: Pre-warm the plate shaker and PCR-thermocycler to 72 °C.

- Prepare PCR-mix with the same N50X for all samples of a plate (note: 96-well plate = almost 2.4 ml of PCR mix!):

<u>24 µl of PCR Mix per one sample:</u>	17 µl	H2O
	5 µl	5x Q5 buffer
	0.5 µl	10 mM dNTP
	1.5 µl	10 uM N50X nextera primer
	0.25 µl	Q5 polymerase

- Remove WASH buffer from the tagmented-DNA-on-beads and place the plate/tubes on ice.

- **Pipette 24 µl of PCR Mix into each well.**

- With a multichannel pipette, **pipette 1.5 µl of 10 uM N70(1-96) nextera primer** (each sample gets different N70X barcoded primer from your barcoded 96-well-N70X-primer plate)

=> PCR volume = 25.5 ul

Mix the plate in plate-mixer preheated to 72°C or by inverting the tubes, and place in a preheated PCR-cycler.

Tip: put the plate on magnetic stand for 2 seconds and then straight into the plate-shaker if the beads have settled already and are difficult to re-suspend just by shaking in plate-shaker.

Important: The 1st step, 10 min at 72°C is absolutely crucial so that the polymerase fills-in the holes made by Tn5. Never heat-denature tagmented samples before this step!

1. 72°C....10 minutes (this fills in 9bp holes made by transposase) (note: **NEVER use HOT-start polymerase**)
2. 98°C....30 sec

4. 98°C... 15 sec

5. 65°C... 20 sec

6. 72°C... 60 sec

7. Go to 4 8x (= 9 cycles in total)

Note: At the end of 10 min 72 °C incubation, pause the cycler and mix the plate in plate-mixer one more time, then un-pause the cycler and continue the protocol.

Spin the plate after PCR to collect all the liquid.

After PCR run 5 µl of several samples on 1.5% agarose gel to see if tagmentation/PCR worked and the size-distribution of the libraries, then add 5 µl of 10mM Tris Buffer back into these samples to make the volume 25 µl again.

Clean-up and size select each sample (see next page) with magnetic beads to get 350-600 bp library -> then pool 5 or 10 ng of each library into the Final Library Pool. Then do one last clean-up with 1:1 ratio of magnetic beads to remove any residual primers and if available run on BioAnalyzer.

Adjust some of the final library pool to 0.8 ng/µl (2.5 nM) to submit at least 15 µl for sequencing (HiSeq300).

SIZE selection of libraries after PCR

1. Add 30 µl of H2O or 10 mM Tris to each of the 25 µl PCR sample to make the volume 55 µl.

Prepare two 96 well plates: **Plate1**. Add 25 µl of size-selection magnetic beads in each well, short spin the plate
Plate2. Add 12 µl of size-selection magnetic beads in each well, short spin the plate

2. Transfer 50 µl of each PCR sample into **Plate1** and mix by pipetting up/down 5-10 times.

- incubate 5 min at RT (*here ration of Beads : Sample = 0.5 => beads bind only DNA above approx. 600+ bp*)

- put on magnet for 3-5 minutes.

3. Transfer the SUPERNATANT (75 µl) into the **Plate2** and mix by pipetting up/down 5-10 times.

- incubate 10-15 min at RT (*here the ration of Beads : Sample increased to 0.74*) => binds everything from approx. 350-600 bp)

-Put on magnet, discard supernatant and wash the beads 2-times with 150 µl of 80% EtOH (not disturbing the beads).

-Remove all the EtOH by pipetting out liquid 2-5 time with the same pipet tip.

-Let stand to evaporate any left-over EtOH (approx.. 2-5 minutes, do not over-dry the beads)

-Take the plate away from the magnet and add 40 µl of 10mM Tris buffer, pH=8.0 to the beads.

-Mix on plate shaker or until the beads re-suspend evenly and incubate 5 min at RT to elute the library DNA.

-Shortly spin the plate down in a centrifuge and put the plate on magnet; let stand 3-5 minutes, and then transfer the supernatant (DNA library) to your final plate.

-Measure concentration of library in 1 or 2 µl of each sample with PicoGreen (Quant-iT PicoGreen dsDNA Assay Kit).

-Pool together 5 ng of each library sample into one 1.5ml eppi to make a "Final Library Pool".

-Check the Final Pool's volume, and do a final clean-up/concentrating of the Final Pool with 1:1 (Pool:Beads) ratio.

If the Final Pool is, let's say, 500 μ l and it is expected to have 480 ng of library DNA (e.g. from 96 samples, 5 ng each), then, to this 500 μ l Final Pool, I would add: 200 μ l of Magnetic-beads + 300ul of Magnetic-Bead-Buffer.

- Vortex and let stand at RT for 10-15 minutes, short spin down and put on the magnet for 3-5 min.
- Remove supernatant; wash the beads on magnet twice with 80% EtOH,
- Evaporate the left-over EtOH for 4 minutes,
- Elute the final library to approximately 1.5-2.5 ng/ μ l expected concentration (in this case, if expecting 480 ng of library, I would use 200 μ l of 10mM Tris buffer, to make a concentrated STOCK-Library pool).
- Adjust the whole pool, or a smaller part (but at least 15 μ l), of the STOCK-Library pool to **0.8 ng/ μ l**, which equals to approximately **2.5nM Library-pool**, required for the sequencing on **HiSeq3000/NovaSeq**.

Buffers / Recipes

DNA extraction/Size-selection-magnetic beads

STOCK BEADS: 15ml, SpeedBeads™ (GE65152105050250 SIGMA) from Sigma /or/ GElifesciences Sera-Mag SpeedBead (65152105050250)

BEAD BUFFER: for 50 ml of Bead buffer, dissolve in 35 ml of mol.biol. grade H₂O:

9 g	PEG8000
7.3 g	NaCl
500 μ l	1M Tris-HCl, pH=8.0
100 μ l	0.5M EDTA, pH=8.0

→ dissolve and add mol.biol. grade H₂O to 50 ml

- on magnetic stand, add 1 ml of the Stock Beads in an 1.5 ml Eppendorf tube
 - after 1 min remove the solution and wash beads with 1 ml of 10mM Tris-HCl, pH=8
 - on magnet, remove Tris buffer and use some of the Bead Buffer to transfer and re-suspend the beads into the rest 50 ml of Bead Buffer.
 - Beads are ready to use, store at 4C, use at RT.
- 1 : 1 ratio (50 μ l sample + 50 μ l beads) => binds all above 200 bp
1 : 0.54 ration (50 μ l sample + 27 μ l beads) => binds all above 500 bp
1 : 0.76 ratio (50 μ l sample + 38 μ l beads) => binds all above 250 pb

5xTn5 buffer = 5x TAPS-Mg-DMF (contains: 50mM TAPS, 25mM MgCl₂, 50% DMF),

DMF = N,N-Dimethylformamide

TAPS = T5130 SIGMA

1M MgCl₂ = e.g. SERVA , 1 M Magnesium chloride, 38772.01

1M TAPS STOCK (50 ml):

1M TAPS in H₂O, adjust to pH:8.5 with 10M NaOH (for 25 ml of 1M TAPS you will need ca 5 ml of 3M NaOH)
⇒ 0.2/0.45 μ m Filter sterilize the 1M TAPS

To prepare 10 ml of 5x Tn5 buffer:

Mix on ice in same order:

4.25 ml H₂O + 500 μ l 1M TAPS (pH=8.5) + 250 μ l 1M MgCl₂ + 5 ml DMF

Note: Buffer becomes warm after adding DMF (keep in ice)

Aliquot in 1.5 ml tubes, store at 4C (for short time, 1-2 months), for longer storage at -20C (buffer doesn't freeze)

Strept. Binding buffer

0.6 M NaCl
10 mM Tris, pH=8.0
0.5 mM EDTA
0.1 % Triton X-100

Dialysis / Tn-bead storage buffer

50 mM HEPES-KOH pH 7.2
0.2 M NaCl
0.2 mM EDTA
2 mM DTT
0.2 % Triton X-100
20 % glycerol

WASH buffer:

10 mM Tris,pH=8,
30 mM NaCl,
0.1 % Triton X-100

“0.6% SDS-Wash”

WASH buffer + 0.6% SDS (from stock 20% SDS in H₂O)

Barcoding oligos

on sample sheet

	<i>If possible: N50X oligos order with PAGE purification and bigger scale especially if planning to pool more than 96 samples per lane into sequencing...</i>	<u>HiSeq3000</u>	<u>NovaSeq</u>
Nextera_N501	AATGATACGGCGACCACCGAGATCTACAC <u>TAGATCGC</u> TCGTCGGCAGCGTC	GCGATCTA	TAGATCGC
Nextera_N502	AATGATACGGCGACCACCGAGATCTACAC <u>CTCTCTAT</u> TCGTCGGCAGCGTC	ATAGAGAG	CTCTCTAT
Nextera_N503	AATGATACGGCGACCACCGAGATCTACAC <u>TATCCTCT</u> TCGTCGGCAGCGTC	AGAGGATA	TATCCTCT
Nextera_N504	AATGATACGGCGACCACCGAGATCTACAC <u>AGAGTAGA</u> TCGTCGGCAGCGTC	TCTACTCT	AGAGTAGA
Nextera_N505	AATGATACGGCGACCACCGAGATCTACAC <u>GTAAGGAG</u> TCGTCGGCAGCGTC	CTCCTTAC	GTAAGGAG
Nextera_N506	AATGATACGGCGACCACCGAGATCTACAC <u>ACTGCATA</u> TCGTCGGCAGCGTC	TATGCACT	ACTGCATA
Nextera_N507	AATGATACGGCGACCACCGAGATCTACAC <u>AAGGAGTA</u> TCGTCGGCAGCGTC	TACTCCTT	AAGGAGTA
Nextera_N508	AATGATACGGCGACCACCGAGATCTACAC <u>CTAAGCCT</u> TCGTCGGCAGCGTC	AGGCTTAG	CTAAGCCT

		barcode in primer		barcode on sample sheet
N701	1	TCGCCTTA	CAAGCAGAAGACGGCATAACGAGAT <u>TCGCCTTA</u> GTCTCGTGGGCTCGG	TAAGGCGA
N702	2	CTAGTACG	CAAGCAGAAGACGGCATAACGAGAT <u>CTAGTACG</u> GTCTCGTGGGCTCGG	CGTACTAG
N703	3	TTCTGCCT	CAAGCAGAAGACGGCATAACGAGAT <u>TTCTGCCT</u> GTCTCGTGGGCTCGG	AGGCAGAA
N704	4	GCTCAGGA	CAAGCAGAAGACGGCATAACGAGAT <u>GCTCAGGA</u> GTCTCGTGGGCTCGG	TCCTGAGC
N705	5	AGGAGTCC	CAAGCAGAAGACGGCATAACGAGAT <u>AGGAGTCC</u> GTCTCGTGGGCTCGG	GGACTCCT
N706	6	CATGCCTA	CAAGCAGAAGACGGCATAACGAGAT <u>CATGCCTA</u> GTCTCGTGGGCTCGG	TAGGCATG
N707	7	GTAGAGAG	CAAGCAGAAGACGGCATAACGAGATGTAGAGAGGTCTCGTGGGCTCGG	CTCTCTAC
N708	8	CCTCTCTG	CAAGCAGAAGACGGCATAACGAGATCCTCTCTGGTCTCGTGGGCTCGG	CAGAGAGG
N709	9	AGCGTAGC	CAAGCAGAAGACGGCATAACGAGATAGCGTAGCGTCTCGTGGGCTCGG	GTACGCT
N710	10	CAGCCTCG	CAAGCAGAAGACGGCATAACGAGATCAGCCTCGGTCTCGTGGGCTCGG	CGAGGCTG
N711	11	TGCCTCTT	CAAGCAGAAGACGGCATAACGAGATTGCCTCTTGTCTCGTGGGCTCGG	AAGAGGCA
N712	12	TCCTCTAC	CAAGCAGAAGACGGCATAACGAGATTCCTCTACGTCTCGTGGGCTCGG	GTAGAGGA
N714	13	TCATAGGC	CAAGCAGAAGACGGCATAACGAGATTCATAGGCGTCTCGTGGGCTCGG	GCCTATGA
N715	14	CCTGAGAT	CAAGCAGAAGACGGCATAACGAGATCCTGAGATGTCTCGTGGGCTCGG	ATCTCAGG
N716	15	TAGCGAGT	CAAGCAGAAGACGGCATAACGAGATTAGCGAGTGTCTCGTGGGCTCGG	ACTCGCTA
N718	16	GTAGCTCC	CAAGCAGAAGACGGCATAACGAGATGTAGCTCCGTCTCGTGGGCTCGG	GGAGCTAC
N719	17	CATCAGGC	CAAGCAGAAGACGGCATAACGAGATCATCAGGCGTCTCGTGGGCTCGG	CGCTAGTG
N720	18	AGGCTCCG	CAAGCAGAAGACGGCATAACGAGATAGGCTCCGGTCTCGTGGGCTCGG	CGGAGCCT
N721	19	GCAGCGTA	CAAGCAGAAGACGGCATAACGAGATGCAGCGTAGTCTCGTGGGCTCGG	TACGCTGC
N722	20	CTGCGCAT	CAAGCAGAAGACGGCATAACGAGATCTGCGCATGTCTCGTGGGCTCGG	ATGCGCAG
N723	21	GAGCGCTA	CAAGCAGAAGACGGCATAACGAGATGAGCGCTAGTCTCGTGGGCTCGG	TAGCGCTC
N724	22	CGCTCAGT	CAAGCAGAAGACGGCATAACGAGATCGCTCAGTGTCTCGTGGGCTCGG	ACTGAGCG
N726	23	TACAGGAT	CAAGCAGAAGACGGCATAACGAGAT <u>TACAGGAT</u> GTCTCGTGGGCTCGG	ATCCTGTA
N727	24	ACTGATCG	CAAGCAGAAGACGGCATAACGAGATACTGATCGGTCTCGTGGGCTCGG	CGATCAGT
N728	25	TAGCTGCA	CAAGCAGAAGACGGCATAACGAGATTAGCTGCAGTCTCGTGGGCTCGG	TGCAGCTA
N729	26	GACGTGCA	CAAGCAGAAGACGGCATAACGAGATGACGTGCGTCTCGTGGGCTCGG	TCGACGTC
N730	27	TGTGGTTG	CAAGCAGAAGACGGCATAACGAGATTGTGGTTGGTCTCGTGGGCTCGG	CAACCACA
N731	28	TAGTCTTG	CAAGCAGAAGACGGCATAACGAGATTAGTCTTGGTCTCGTGGGCTCGG	CAAGACTA
N732	29	TTCCATTG	CAAGCAGAAGACGGCATAACGAGATTTCCATTGGTCTCGTGGGCTCGG	CAATGGAA
N733	30	TCGAAGTG	CAAGCAGAAGACGGCATAACGAGATTCGAAGTGGTCTCGTGGGCTCGG	CACCTCGA
N734	31	TAACGCTG	CAAGCAGAAGACGGCATAACGAGATTAACGCTGGTCTCGTGGGCTCGG	CAGCGTTA
N735	32	TTGGTATG	CAAGCAGAAGACGGCATAACGAGATTTGGTATGGTCTCGTGGGCTCGG	CATACCAA
N736	33	TGAAGTGG	CAAGCAGAAGACGGCATAACGAGATTGAACTGGTCTCGTGGGCTCGG	CCAGTTCA
N737	34	TAGTTCCG	CAAGCAGAAGACGGCATAACGAGATTAGTTCGGGTCTCGTGGGCTCGG	CCGAACATA
N738	35	TCTCACGG	CAAGCAGAAGACGGCATAACGAGATTCACGGGTCTCGTGGGCTCGG	CCGTGAGA
N739	36	TCATGAGG	CAAGCAGAAGACGGCATAACGAGATTCATGAGGTCTCGTGGGCTCGG	CCTCATGA
N740	37	TAAGTTCG	CAAGCAGAAGACGGCATAACGAGATTAAGTTCGGTCTCGTGGGCTCGG	CGAACTTA
N741	38	TCCAGTCC	CAAGCAGAAGACGGCATAACGAGATTCAGTCCGGTCTCGTGGGCTCGG	CGACTGGA
N742	39	TGTATGCG	CAAGCAGAAGACGGCATAACGAGATTGTATGCGGTCTCGTGGGCTCGG	CGCATACA
N743	40	TCATTGAG	CAAGCAGAAGACGGCATAACGAGATTCATTGAGGTCTCGTGGGCTCGG	CTCAATGA
N744	40	TGACTCAG	CAAGCAGAAGACGGCATAACGAGATTGACTCAGGTCTCGTGGGCTCGG	CTGAGTCA
N745	42	TATGCCAG	CAAGCAGAAGACGGCATAACGAGATTTATGCCAGGTCTCGTGGGCTCGG	CTGGCATA
N746	43	TCAGATTC	CAAGCAGAAGACGGCATAACGAGATTCAGATTCGATTCGTTCTCGTGGGCTCGG	GAATCTGA
N747	44	TACTAGTC	CAAGCAGAAGACGGCATAACGAGATTAAGTTCGTTCTCGTGGGCTCGG	ACTAGTA
N748	45	TTCAGCTC	CAAGCAGAAGACGGCATAACGAGATTTTCAGCTCGTCTCGTGGGCTCGG	GAGCTGAA
N749	46	TGTCTATC	CAAGCAGAAGACGGCATAACGAGATTGTCTATCGTCTCGTGGGCTCGG	GATAGACA
N750	47	TATGTGGC	CAAGCAGAAGACGGCATAACGAGATTAATGTGGGTCTCGTGGGCTCGG	GCCACATA
N751	48	TTACTCGC	CAAGCAGAAGACGGCATAACGAGATTTACTCGGTCTCGTGGGCTCGG	GCGAGTAA
N752	49	TCGTTAGC	CAAGCAGAAGACGGCATAACGAGATTCGTTAGCTCTCGTGGGCTCGG	GCTAACGA

N753	50	TACGGAGC	CAAGCAGAAGACGGCATAACGAGATTACGGAGCGTCTCGTGGGCTCGG	GCTCCGTA
N754	51	TGTTCTCC	CAAGCAGAAGACGGCATAACGAGATTGTTCTCCGTCTCGTGGGCTCGG	GGAGAACA
N755	52	TTCGCACC	CAAGCAGAAGACGGCATAACGAGATTTTCGCACCGTCTCGTGGGCTCGG	GGTGCGAA
N756	53	TTGCGTAC	CAAGCAGAAGACGGCATAACGAGATTTGCGTACGTCTCGTGGGCTCGG	GTACGCAA
N757	54	TCTACGAC	CAAGCAGAAGACGGCATAACGAGATTTCTACGACGTCTCGTGGGCTCGG	GTCTGTAGA
N758	55	TGACAGAC	CAAGCAGAAGACGGCATAACGAGATTTGACAGACGTCTCGTGGGCTCGG	GTCTGTCA
N759	55	TAGAAAAC	CAAGCAGAAGACGGCATAACGAGATTTAGAAAACGTCTCGTGGGCTCGG	GTTTTCTA
N760	57	TCATCCTA	CAAGCAGAAGACGGCATAACGAGATTCATCCTAGTCTCGTGGGCTCGG	TAGGATGA
N761	58	TGCTGATA	CAAGCAGAAGACGGCATAACGAGATTGCTGATAGTCTCGTGGGCTCGG	TATCAGCA
N762	59	TAGACGGA	CAAGCAGAAGACGGCATAACGAGATTAGACGGAGTCTCGTGGGCTCGG	TCCGTCTA
N763	60	TGTGAAGA	CAAGCAGAAGACGGCATAACGAGATTTGTGAAGAGTCTCGTGGGCTCGG	TCTTCACA
N764	61	TCTCTTCA	CAAGCAGAAGACGGCATAACGAGATTTCTTTCAGTCTCGTGGGCTCGG	TGAAGAGA
N765	62	TTGTTCCA	CAAGCAGAAGACGGCATAACGAGATTTTGTCCAGTCTCGTGGGCTCGG	TGGAACAA
N766	63	TGAAGCCA	CAAGCAGAAGACGGCATAACGAGATTTGAAGCCAGTCTCGTGGGCTCGG	TGGCTTCA
N767	64	TACCACCA	CAAGCAGAAGACGGCATAACGAGATTACCACCAGTCTCGTGGGCTCGG	TGGTGGTA
N768	65	TGCGTGAA	CAAGCAGAAGACGGCATAACGAGATTGCGTGAAAGTCTCGTGGGCTCGG	TTACGCA
N769	66	GGTGAGTT	CAAGCAGAAGACGGCATAACGAGATGGTGAGTTGTCTCGTGGGCTCGG	AACTCACC
N770	67	GATCTCTT	CAAGCAGAAGACGGCATAACGAGATGATCTCTTGTCTCGTGGGCTCGG	AAGAGATC
N771	68	GTGTCCCT	CAAGCAGAAGACGGCATAACGAGATGTGTCTTGTCTCGTGGGCTCGG	AAGGACAC
N772	69	GACGGATT	CAAGCAGAAGACGGCATAACGAGATGACGGATTGTCTCGTGGGCTCGG	AATCCGTC
N773	70	GCAACATT	CAAGCAGAAGACGGCATAACGAGATGCAACATTGTCTCGTGGGCTCGG	AATGTTGC
N774	71	GGTCGTGT	CAAGCAGAAGACGGCATAACGAGATGGTCGTGTGTCTCGTGGGCTCGG	ACACGACC
N775	72	GAATCTGT	CAAGCAGAAGACGGCATAACGAGATGAATCTGTGTCTCGTGGGCTCGG	ACAGATTC
N776	73	GTACATCT	CAAGCAGAAGACGGCATAACGAGATGTACATCTGTCTCGTGGGCTCGG	AGATGTAC
N777	74	GAGGTGCT	CAAGCAGAAGACGGCATAACGAGATGAGGTGCTGTCTCGTGGGCTCGG	AGCACCTC
N778	75	GCATGGCT	CAAGCAGAAGACGGCATAACGAGATGCATGGCTGTCTCGTGGGCTCGG	AGCCATGC
N779	76	GTTAGCCT	CAAGCAGAAGACGGCATAACGAGATGTTAGCCTGTCTCGTGGGCTCGG	AGGCTAAC
N780	77	GTCGCTAT	CAAGCAGAAGACGGCATAACGAGATGTCGCTATGTCTCGTGGGCTCGG	ATAGCGAC
N781	78	GGAATGAT	CAAGCAGAAGACGGCATAACGAGATGGAATGATGTCTCGTGGGCTCGG	ATCATTCC
N782	79	GAGCCAAT	CAAGCAGAAGACGGCATAACGAGATGAGCCAATGTCTCGTGGGCTCGG	ATTGGCTC
N783	80	GCTCCTTG	CAAGCAGAAGACGGCATAACGAGATGCTCCTTGGTCTCGTGGGCTCGG	CAAGGAGC
N784	81	GTAAGGTG	CAAGCAGAAGACGGCATAACGAGATGTAAGGTGGTCTCGTGGGCTCGG	CACCTTAC
N785	82	GAGGATGG	CAAGCAGAAGACGGCATAACGAGATGAGGATGGGTCTCGTGGGCTCGG	CCATCCTC
N786	83	GTTGTCCG	CAAGCAGAAGACGGCATAACGAGATGTTGTCCGGTCTCGTGGGCTCGG	CCGACAAC
N787	84	GGATTAGG	CAAGCAGAAGACGGCATAACGAGATGGATTAGGGTCTCGTGGGCTCGG	CCTAATCC
N788	85	GATAGAGG	CAAGCAGAAGACGGCATAACGAGATGATAGAGGGTCTCGTGGGCTCGG	CCTCTATC
N789	86	GTGTGTCC	CAAGCAGAAGACGGCATAACGAGATGTGTGTCCGGTCTCGTGGGCTCGG	CGACACAC
N790	87	GCAATCCG	CAAGCAGAAGACGGCATAACGAGATGCAATCCGGTCTCGTGGGCTCGG	CGGATTGC
N791	88	GACCTTAG	CAAGCAGAAGACGGCATAACGAGATGACCTTAGGTCTCGTGGGCTCGG	CTAAGGTC
N792	89	GCCTGTTC	CAAGCAGAAGACGGCATAACGAGATGCCTGTTCGTCTCGTGGGCTCGG	GAACAGGC
N793	90	GCACTGTC	CAAGCAGAAGACGGCATAACGAGATGCACTGTCTCGTGGGCTCGG	GACAGTGC
N794	91	GCTAACTC	CAAGCAGAAGACGGCATAACGAGATGCTAACTCGTCTCGTGGGCTCGG	GAGTTAGC
N795	92	GATTCATC	CAAGCAGAAGACGGCATAACGAGATGATTCATCGTCTCGTGGGCTCGG	GATGAATC
N796	93	GTCTTGCC	CAAGCAGAAGACGGCATAACGAGATGTCTTGGCTCTCGTGGGCTCGG	GCCAAGAC
N797	94	TGCGATCT	CAAGCAGAAGACGGCATAACGAGATTTGCGATCTGTCTCGTGGGCTCGG	AGATCGCA
N798	95	TTCCTGCT	CAAGCAGAAGACGGCATAACGAGATTTCTGCTGTCTCGTGGGCTCGG	AGCAGGAA
N799	96	TAGACACT	CAAGCAGAAGACGGCATAACGAGATTAGACACTGTCTCGTGGGCTCGG	AGTGTCTA

NOTES:

N50X primers with PAGE purification can prevent data loss if multiplexing more than 96 samples, with several different N50X primers. In that case, you have to also use the N50X barcode to demultiplex data, not just the N70X barcode. PAGE purification keeps the errors in the barcode at minimum and thus helps with data recovery. It cost more and one needs to order bigger scale, but it is worth it. Of course, ideal would be to order not just N50X primers with PAGE purification, but also the N70X primers, but that would cost a lot.

Normally, if you have less than 96 samples in a sequencing pool, you can use the cheap, NOT-PAGE-purified-N50X, because you can ignore N50X barcode during demultiplexing (all samples have the same N50X) and thus you can be de-multiplex only using N70X barcode and recover more data (e.g. you will not lose any reads due to errors in N50X barcode, which you would, if you demultiplex using both N70X and N50X barcodes).